

Product Datasheet

HTRA2 (Phospho-S142) antibody (orb571001)

Description	HTRA2 (Phospho-S142) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human HtrA2 around the phosphorylation site of S142.
Target	HTRA2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O43464
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:40000
Expiration Date	12 months from date of receipt.