

Product Datasheet

EEF2K (Phospho-S359) antibody (orb571015)

Description	EEF2K (Phospho-S359) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human eEF2K around the phosphorylation site of S359.
Target	EEF2K
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	000418
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:20000
Expiration Date	12 months from date of receipt.