

Product Datasheet

HRH1 (Phospho-S398) antibody (orb571031)

Description	HRH1 (Phospho-S398) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Synthesized peptide derived from Human Histamine H1 Receptor around the phosphorylation site of S398.
Target	HRH1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35367
Dilution Range	IF: 1:200-1:1000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.