

Product Datasheet

IGF1R (Phospho-Y1346) antibody (orb571041)

Description	IGF1R (Phospho-Y1346) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human IGF-IR around the phosphorylation site of Y1346.
Target	IGF1R
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08069
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:40000
Expiration Date	12 months from date of receipt.