

Product Datasheet

FAS (Phospho-Y291) antibody (orb571044)

Description	FAS (Phospho-Y291) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Synthesized peptide derived from Human FAS around the phosphorylation site of Y291.
Target	FAS
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25445
Dilution Range	IF: 1:200-1:1000, ELISA: 1:5000
Expiration Date	12 months from date of receipt.