

## Product Datasheet

### MARK1 (Phospho-T215) antibody (orb571090)

|                            |                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------|
| <b>Description</b>         | MARK1 (Phospho-T215) antibody                                                                     |
| <b>Species/Host</b>        | Rabbit                                                                                            |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                      |
| <b>Tested Applications</b> | ELISA, IHC                                                                                        |
| <b>Immunogen</b>           | Synthesized peptide derived from Human MARK1/2/3/4 around the phosphorylation site of T215.       |
| <b>Target</b>              | MARK1/MARK2/MARK3/MARK4                                                                           |
| <b>Preservatives</b>       | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.                      |
| <b>Form/Appearance</b>     | Liquid                                                                                            |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                     |
| <b>Note</b>                | For research use only                                                                             |
| <b>Isotype</b>             | IgG                                                                                               |
| <b>Clonality</b>           | Polyclonal                                                                                        |
| <b>Uniprot ID</b>          | <a href="#">P27448</a> , <a href="#">Q9P0L2</a> , <a href="#">Q7KZI7</a> , <a href="#">Q96L34</a> |
| <b>Dilution Range</b>      | IHC-P: 1:100-1:300, ELISA: 1:5000                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                   |