

Product Datasheet

MYC (Phospho-S373) antibody (orb571112)

Description	MYC (Phospho-S373) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP
Immunogen	Synthesized peptide derived from Human c-Myc around the phosphorylation site of S373.
Target	MYC
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01106
Dilution Range	IHC-P: 1:100-1:300, IP: 2-5 µg: mg lysate, ELISA: 1:20000
Expiration Date	12 months from date of receipt.