

Product Datasheet

NOS2 (Phospho-Y151) antibody (orb571138)

Description	NOS2 (Phospho-Y151) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human NOS2 around the phosphorylation site of Y151.
Target	NOS2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35228
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:10000
Expiration Date	12 months from date of receipt.