

Product Datasheet

NTRK1 (Phospho-Y496) antibody (orb571148)

Description	NTRK1 (Phospho-Y496) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Synthesized peptide derived from Human Trk A around the phosphorylation site of Y496.
Target	NTRK1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04629
Dilution Range	IF: 1:200-1:1000, ELISA: 1:5000
Expiration Date	12 months from date of receipt.