

Product Datasheet

RBL1 (Phospho-T369) antibody (orb571248)

Description	RBL1 (Phospho-T369) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human p107 around the phosphorylation site of T369.
Target	RBL1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P28749
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:5000
Expiration Date	12 months from date of receipt.