

Product Datasheet

RELA (Phospho-T435) antibody (orb571251)

Description	RELA (Phospho-T435) antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Synthesized peptide derived from Human NFκB-p65 around the phosphorylation site of T435.
Target	RELA
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q04206
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:100-1:300, IP: 2-5 µg: mg lysate, ELISA: 1:20000
Expiration Date	12 months from date of receipt.