


Product Datasheet

ABPP (phospho-Thr668) antibody (orb5713)

Description	ABPP (phospho-Thr668) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human APP around the phosphorylation site of Thr668 (AV(p-T)PE)
Target	phospho-Amyloid Precursor Protein (Thr668)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	3.4-9/83 kDa
Uniprot ID	P05067
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Human SH-SY5Y cell; L...