

Product Datasheet

TCF3 (Phospho-T355) antibody (orb571327)

Description	TCF3 (Phospho-T355) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human E2A around the phosphorylation site of T355.
Target	TCF3
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15923
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:5000
Expiration Date	12 months from date of receipt.