

Product Datasheet

TGFBR2 (Phospho-S225) antibody (orb571333)

Description	TGFBR2 (Phospho-S225) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human TGFβ RII around the phosphorylation site of S225.
Target	TGFBR2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P37173
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:10000
Expiration Date	12 months from date of receipt.