

Product Datasheet

VHL (Phospho-S68) antibody (orb571362)

Description	VHL (Phospho-S68) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human VHL around the phosphorylation site of S68.
Target	VHL
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P40337
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:5000
Expiration Date	12 months from date of receipt.