

Product Datasheet

CASP2 (Phospho-S157) antibody (orb571383)

Description	CASP2 (Phospho-S157) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human Caspase-2 around the phosphorylation site of S157.
Target	CASP2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P42575
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:10000
Expiration Date	12 months from date of receipt.