

Product Datasheet

ITCH (Phospho-Y420) antibody (orb571386)

Description	ITCH (Phospho-Y420) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human AIP4 around the phosphorylation site of Y420.
Target	ITCH
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96J02
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:20000
Expiration Date	12 months from date of receipt.