

Product Datasheet

AR (Phospho-S94) antibody (orb571471)

Description	AR (Phospho-S94) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Synthesized peptide derived from Human AR around the phosphorylation site of S94.
Target	AR
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10275
Dilution Range	IF: 1:200-1:1000, ELISA: 1:5000
Expiration Date	12 months from date of receipt.