

Product Datasheet

GSK3B (Phospho-S9) antibody (orb571541)

Description	GSK3B (Phospho-S9) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Synthesized peptide derived from Human GSK3 β around the phosphorylation site of S9.
Target	GSK3B
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49841
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:100-1:300, IP: 2-5 μ g: mg lysate, ELISA: 1:5000
Expiration Date	12 months from date of receipt.