

Product Datasheet

MAP2K2 (Phospho-T394) antibody (orb571559)

Description	MAP2K2 (Phospho-T394) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Synthesized peptide derived from Human MEK-2 around the phosphorylation site of T394.
Target	MAP2K2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P36507
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:100-1:300, IP: 2-5 µg: mg lysate, ELISA: 1:10000
Expiration Date	12 months from date of receipt.