

Product Datasheet

MDM2 (Phospho-Tyr394) antibody (orb571651)

Description	MDM2 (Phospho-Tyr394) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide from Human protein at AA range: 330-410
Target	MDM2
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q00987
Dilution Range	WB: 1:500-2000, ELISA: 1:10000-20000
Expiration Date	12 months from date of receipt.