

Product Datasheet

TH (Phospho-Ser40) antibody (orb571683)

Description	TH (Phospho-Ser40) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Peptide sequence around phosphorylation site of serine 40 (R-Q-S(p)-L-I) derived from Human Tyrosine Hydroxylase (TH).
Target	TH
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07101
Dilution Range	IF: 1:100-1:200
Expiration Date	12 months from date of receipt.