

Product Datasheet

EIF2S2 (Phospho-Ser2) antibody (orb571710)

Description	EIF2S2 (Phospho-Ser2) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of serine 2 (M-S(p)-G-D-E) derived from Human eIF2?.
Target	EIF2S2
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20042
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.