

## Product Datasheet

### PRKCZ (Phospho-Thr410) antibody (orb571734)

|                            |                                                                                                                                     |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | PRKCZ (Phospho-Thr410) antibody                                                                                                     |
| <b>Species/Host</b>        | Rabbit                                                                                                                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                        |
| <b>Tested Applications</b> | ELISA, IHC                                                                                                                          |
| <b>Immunogen</b>           | Peptide sequence around phosphorylation site of threonine 410(T-S-T(p)-F-C) derived from Human PKCζ.                                |
| <b>Target</b>              | PRKCZ                                                                                                                               |
| <b>Preservatives</b>       | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| <b>Form/Appearance</b>     | liquid                                                                                                                              |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                       |
| <b>Note</b>                | For research use only                                                                                                               |
| <b>Isotype</b>             | IgG                                                                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                                                                          |
| <b>Uniprot ID</b>          | <a href="#">Q05513</a>                                                                                                              |
| <b>Dilution Range</b>      | IHC-P: 1:50-1:100                                                                                                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                     |