

Product Datasheet

PTPN1 (Phospho-Ser50) antibody (orb571768)

Description	PTPN1 (Phospho-Ser50) antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of Serine 50(D-V-S(p)-P-F) derived from Human PTP1B.
Target	PTPN1
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P18031
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.