

Product Datasheet

KRT8 (Phospho-Ser432) antibody (orb571771)

Description	KRT8 (Phospho-Ser432) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of Serine 432(L-T-S(p)-P-G) derived from Human Keratin 8.
Target	KRT8
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05787
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.