

Product Datasheet

GATA4 (Phospho-Ser262) antibody (orb571772)

Description	GATA4 (Phospho-Ser262) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of Serine 262(R-L-S(p)-A-S) derived from Human GATA4.
Target	GATA4
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P43694
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.