

Product Datasheet

ATP1A1 (Phospho-Ser23) antibody (orb571781)

Description	ATP1A1 (Phospho-Ser23) antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of Serine 23 (K-K-G(p)-K-K) derived from Rat Na,K-ATPase alpha1.
Target	ATP1A1
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05023
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.