

Product Datasheet

GFAP (Phospho-Ser38) antibody (orb571785)

Description	GFAP (Phospho-Ser38) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of Serine 38(R-L-S(p)-L-A) derived from Human GFAP.
Target	GFAP
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14136
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.