

Product Datasheet

PRKD1 (Phospho-Ser738/742) antibody (orb571803)

Description	PRKD1 (Phospho-Ser738/742) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of Serine738/Serine 742(E-K-S(p)-F-R-R-S(p)-V-V) derived from Human PKD1/2/3/PKC μ .
Target	PRKD1/PRKD2/PRKD3
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15139, O94806, Q9BZL6
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.