

Product Datasheet

CD19 (Phospho-Tyr531) antibody (orb571816)

Description	CD19 (Phospho-Tyr531) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of tyrosine 531 (D-S-Y(p)-E-N) derived from Human CD19.
Target	CD19
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15391
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.