

Product Datasheet

RPS6KA5 (Phospho-Ser212) antibody (orb571835)

Description	RPS6KA5 (Phospho-Ser212) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of serine 212 (A-Y-S(p)-F-C) derived from Human MSK1.
Target	RPS6KA5
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	075582
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.