

Product Datasheet

CDC25B (Phospho-Ser353) antibody (orb571847)

Description	CDC25B (Phospho-Ser353) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of serine 353(R-R-S(p)-V-T) derived from Human CDC25B .
Target	CDC25B
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P30305
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.