

Product Datasheet

PTK2 (Phospho-Tyr576) antibody (orb571854)

Description	PTK2 (Phospho-Tyr576) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of Tyrosine576 (S-T-Y(p)-Y-K) derived from Human FAK.
Target	PTK2
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q05397
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.