

Product Datasheet

GATA2 (Phospho-Ser401) antibody (orb571903)

Description	GATA2 (Phospho-Ser401) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of Serine 401(K-M-S(p)-N-K) derived from Human GATA2 .
Target	GATA2
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P23769
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.