

Product Datasheet

TGFR2 (Phospho-Ser225/250) antibody (orb571907)

Description	Rabbit polyclonal antibody to TGFR2 (Phospho-Ser225/250)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Peptide sequence around phosphorylation site of Serine 225/250(D-R-S(p)-D-I) derived from Human TGF β Receptor II.
Target	TGFR2
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P37173
Dilution Range	IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.