

Product Datasheet

MAP2K7 (Phospho-Thr275) antibody (orb571926)

Description	MAP2K7 (Phospho-Thr275) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of threonine 275(A-K-T(p)-R-S) derived from Human MAP2K7.
Target	MAP2K7
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14733
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.