

Product Datasheet

PKMYT1 (Phospho-Ser83) antibody (orb571927)

Description	PKMYT1 (Phospho-Ser83) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Peptide sequence around phosphorylation site of Serine 83(R-V-S(p)-F-R) derived from Human MYT1.
Target	PKMYT1
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99640
Dilution Range	IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.