

Product Datasheet

TUBA1A (Phospho-Tyr272) antibody (orb571982)

Description	TUBA1A (Phospho-Tyr272) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around phosphorylation site of tyrosine 272 (A-T-Y(p)-A-P) derived from Human TUBA1/3/4.
Target	TUBA1A/TUBA1B/TUBA1C/TUBA3C/TUBA3E
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q71U36
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.