

Product Datasheet

JUNB (Phospho-Ser79) antibody (orb572037)

Description	JUNB (Phospho-Ser79) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Peptide sequence around phosphorylation site of serine 79 (G-A-S(p)-L-K) derived from Human JunB.
Target	JUNB
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P17275
Dilution Range	IHC-P: 1:50-1:100, IF: 1:100-1:200
Expiration Date	12 months from date of receipt.