

Product Datasheet

PLCG2 (Phospho-Tyr753) antibody (orb572097)

Description	PLCG2 (Phospho-Tyr753) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Peptide sequence around phosphorylation site of tyrosine 753 (S-L-Y(p)-D-V) derived from Human PLCg2.
Target	PLCG2
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16885
Dilution Range	WB: 1:500-1:1000, IF: 1:100-1:200
Expiration Date	12 months from date of receipt.