

Product Datasheet

Arfaptin 2 (phospho-S260) antibody (orb5722)

Description	Arfaptin 2 (phospho-S260) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human arfaptin 2 around the phosphorylation site of Ser260 (LE(p-S)AQ)
Target	phospho-arfaptin 2 (Ser260)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38 kDa
Uniprot ID	P53365
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.