

Product Datasheet

MAPK8 (Phospho-Thr183/Tyr185) antibody (orb572314)

Description	MAPK8 (Phospho-Thr183/Tyr185) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Peptide sequence around phosphorylation site of Thr183/Tyr185 (M-M-T(p)-P-Y(p)- V - V) derived from Human JNK1/JNK2/JNK3.
Target	MAPK8/MAPK9/MAPK10
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P53779, P45983, P45984
Dilution Range	WB: 1:500-1:1000, IF: 1:100-1:200
Expiration Date	12 months from date of receipt.