

Product Datasheet

MUC1 (Phospho-Tyr1229) antibody (orb572337)

Description	MUC1 (Phospho-Tyr1229) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Peptide sequence around phosphorylation site of tyrosine 1229 S-P-Y(p)-E-K) derived from Human CD227/MUC1.
Target	MUC1
Preservatives	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15941
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, IF: 1:100-1:200
Expiration Date	12 months from date of receipt.