

Product Datasheet

GSK3A (Phospho-Tyr279/216) antibody (orb572363)

Description	GSK3A (Phospho-Tyr279/216) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Peptide sequence around phosphorylation site of tyrosine 279/216 (V-S-Y(p)-I-C) derived from Human GSK3α/β .
Target	GSK3A/GSK3B
Preservatives	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49840, P49841
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, IF: 1:100-1:200
Expiration Date	12 months from date of receipt.