

Product Datasheet

Mono-Methyl-Histone H1 (Lys25) antibody (orb572388)

Description	Mono-Methyl-Histone H1 (Lys25) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic Peptide
Target	HIST1H1A/HIST1H1C/HIST1H1D/HIST1H1E
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10412, P16402, Q02539, P16403
Dilution Range	WB: 1:500-1:1000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.