

Product Datasheet

Mono-Methyl-Histone H2B (Lys5) antibody (orb572390)

Description	Mono-Methyl-Histone H2B (Lys5) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic Peptide
Target	HIST1H2BA/HIST1H2BB/HIST1H2BC
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P33778, P62807, Q96A08
Dilution Range	WB: 1:500-1:1000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.