

Product Datasheet

Tri-Methyl-Histone H3 (Lys79) antibody (orb572433)

Description	Tri-Methyl-Histone H3 (Lys79) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic Peptide
Target	HIST1H3A/HIST2H3A/H3F3A
Preservatives	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q71DI3 , P84243 , P68431
Dilution Range	WB: 1:500-1:1000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.