

Product Datasheet

Histone H3K79me3 antibody (orb572455)

Description	Histone H3K79me3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IF, IHC, IP, WB
Immunogen	A synthetic methylated peptide corresponding to residues surrounding K79 of Human histone H3
Target	HIST3H3
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16695
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200, IP: 1:50-1:200, ChIP: 1:50-1:200
Expiration Date	12 months from date of receipt.