

Product Datasheet

Histone H3R2me2s antibody (orb572456)

Description	Histone H3R2me2s antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	A synthetic methylated peptide corresponding to residues surrounding R2 of Human histone H3
Target	HIST3H3
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16695
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200, IP: 1:50-1:200
Expiration Date	12 months from date of receipt.